

Product datasheet for **GA117778**

Human **LRRC37A3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	LRRC37A3 CRISPRa kit - CRISPR gene activation of human leucine rich repeat containing 37 member A3
Symbol:	LRRC37A3
UniProt ID:	O60309
Synonyms:	LRRC37; LRRC37A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)