

Product datasheet for **GA117784**

Human TEX45 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TEX45 CRISPRa kit - CRISPR gene activation of human testis expressed 45
Symbol:	TEX45
Components:	GA117784G1, TEX45 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117784G2, TEX45 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117784G3, TEX45 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8NA69</u>
Synonyms:	FLJ35784; FLJ56642
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)