


Product datasheet for **GA117794**

Human ZNF773 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF773 CRISPRa kit - CRISPR gene activation of human zinc finger protein 773
Symbol:	ZNF773
Components:	GA117794G1, ZNF773 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117794G2, ZNF773 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117794G3, ZNF773 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6PK81</u>
Synonyms:	ZNF419B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)