

Product datasheet for **GA117798**

Human C1orf223 (TEX38) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TEX38 CRISPRa kit - CRISPR gene activation of human testis expressed 38
Symbol:	TEX38
Components:	GA117798G1, C1orf223 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117798G2, C1orf223 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117798G3, C1orf223 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6PEX7</u>
Synonyms:	ATPAF1-AS1; C1orf223; THEG4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)