

Product datasheet for **GA117806**

Human Profilin 4 (PFN4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PFN4 CRISPRa kit - CRISPR gene activation of human profilin family member 4
Symbol:	PFN4
UniProt ID:	<u>Q8NHR9</u>
Synonyms:	profilin family, member 4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)