

Product datasheet for **GA117808**

Human ANKRD36 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ANKRD36 CRISPRa kit - CRISPR gene activation of human ankyrin repeat domain 36
Symbol:	ANKRD36
Components:	GA117808G1, ANKRD36 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117808G2, ANKRD36 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117808G3, ANKRD36 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A6QL64
Synonyms:	UNQ2430
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)