

Product datasheet for **GA117828**

Human KCP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCP CRISPRa kit - CRISPR gene activation of human kielin cysteine rich BMP regulator
Symbol:	KCP
Components:	GA117828G1, KCP gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117828G2, KCP gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117828G3, KCP gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZWI8</u>
Synonyms:	CRIM2; KCP1; NET67
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)