

Product datasheet for **GA117833**

Human SWI5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SWI5 CRISPRa kit - CRISPR gene activation of human SWI5 homologous recombination repair protein
Symbol:	SWI5
Components:	GA117833G1, SWI5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117833G2, SWI5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117833G3, SWI5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q1ZZU3
Synonyms:	C9orf119; SAE3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)