

Product datasheet for **GA117838**

Human **LRRC10** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	LRRC10 CRISPRa kit - CRISPR gene activation of human leucine rich repeat containing 10
Symbol:	LRRC10
Components:	GA117838G1, LRRC10 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117838G2, LRRC10 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117838G3, LRRC10 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5BKY1</u>
Synonyms:	HRLRRP; LRRC10A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)