

Product datasheet for **GA117856**

Human KCTD4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCTD4 CRISPRa kit - CRISPR gene activation of human potassium channel tetramerization domain containing 4
Symbol:	KCTD4
Components:	GA117856G1, KCTD4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117856G2, KCTD4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117856G3, KCTD4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8WVF5</u>
Synonyms:	bA321C24.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)