

Product datasheet for **GA117873**

Human ZNF307 (ZKSCAN4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZKSCAN4 CRISPRa kit - CRISPR gene activation of human zinc finger with KRAB and SCAN domains 4
Symbol:	ZKSCAN4
UniProt ID:	Q969J2
Synonyms:	P1P373C6; ZNF307; ZNF427; ZSCAN36
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)