

Product datasheet for **GA117876**

Human C6orf174 (SOGA3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SOGA3 CRISPRa kit - CRISPR gene activation of human SOGA family member 3
Symbol:	SOGA3
Components:	GA117876G1, C6orf174 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117876G2, C6orf174 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117876G3, C6orf174 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5TF21</u>
Synonyms:	C6orf174; dj403A15.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)