

Product datasheet for **GA117894**

Human C10orf140 (SKIDA1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SKIDA1 CRISPRa kit - CRISPR gene activation of human SKI/DACH domain containing 1
Symbol:	SKIDA1
Components:	GA117894G1, C10orf140 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117894G2, C10orf140 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117894G3, C10orf140 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q1XH10</u>
Synonyms:	C10orf140; DLN-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)