

Product datasheet for **GA117902**

Human TEX36 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TEX36 CRISPRa kit - CRISPR gene activation of human testis expressed 36
Symbol:	TEX36
Components:	GA117902G1, TEX36 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117902G2, TEX36 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117902G3, TEX36 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5VZQ5</u>
Synonyms:	bA383C5.1; C10orf122
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)