


## Product datasheet for **GA117910**

### Human Speedy protein C (SPDYC) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	SPDYC CRISPRa kit - CRISPR gene activation of human speedy/RINGO cell cycle regulator family member C
Symbol:	SPDYC
UniProt ID:	<a href="#">Q5MJ68</a>
Synonyms:	Ringo2; RINGOC
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)