

## Product datasheet for **GA117919**

### Human TMEM233 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TMEM233 CRISPRa kit - CRISPR gene activation of human transmembrane protein 233
Symbol:	TMEM233
UniProt ID:	<a href="#">B4DJY2</a>
Synonyms:	DSPB2; IFITMD2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)