


## Product datasheet for **GA117923**

### Human NHLRC3 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | NHLRC3 CRISPRa kit - CRISPR gene activation of human NHL repeat containing 3                                                                                                                                                                                                                                  |
| Symbol:       | NHLRC3                                                                                                                                                                                                                                                                                                        |
| Components:   | GA117923G1, NHLRC3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA117923G2, NHLRC3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA117923G3, NHLRC3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q5JS37</a>                                                                                                                                                                                                                                                                                        |
| Synonyms:     | DKFZp313M1221; DKFZp686E1140                                                                                                                                                                                                                                                                                  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)