

Product datasheet for **GA117927**

Human TMEM179 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM179 CRISPRa kit – CRISPR gene activation of human transmembrane protein 179
Symbol:	TMEM179
Components:	GA117927G1, TMEM179 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117927G2, TMEM179 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117927G3, TMEM179 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZVK1</u>
Synonyms:	C14orf90; TMEM179A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

