


Product datasheet for **GA117930**

Human TNFAIP8L3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNFAIP8L3 CRISPRa kit - CRISPR gene activation of human TNF alpha induced protein 8 like 3
Symbol:	TNFAIP8L3
Components:	GA117930G1, TNFAIP8L3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117930G2, TNFAIP8L3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117930G3, TNFAIP8L3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q5GJ75
Synonyms:	TIPE3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)