

Product datasheet for **GA117935**

Human C16orf87 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C16orf87 CRISPRa kit - CRISPR gene activation of human chromosome 16 open reading frame 87
Symbol:	C16orf87
Components:	GA117935G1, C16orf87 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117935G2, C16orf87 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117935G3, C16orf87 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6PH81</u>
Synonyms:	MGC62018
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)