

Product datasheet for **GA117944**

Human SPDYE4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SPDYE4 CRISPRa kit - CRISPR gene activation of human speedy/RINGO cell cycle regulator family member E4
Symbol:	SPDYE4
Components:	GA117944G1, SPDYE4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117944G2, SPDYE4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117944G3, SPDYE4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A6NLX3</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)