

Product datasheet for **GA117945**

Human TMEM220 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM220 CRISPRa kit - CRISPR gene activation of human transmembrane protein 220
Symbol:	TMEM220
Components:	GA117945G1, TMEM220 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117945G2, TMEM220 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117945G3, TMEM220 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6QAJ8</u>
Synonyms:	A730055C05Rik; MGC130169; MGC130170; OTTMUSP00000006112; transmembrane protein 220
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)