

Product datasheet for **GA117947**

Human C17orf76 (LRRC75A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LRRC75A CRISPRa kit - CRISPR gene activation of human leucine rich repeat containing 75A
Symbol:	LRRC75A
UniProt ID:	<u>Q8NAA5</u>
Synonyms:	C17orf76; FAM211A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)