

Product datasheet for **GA117948**

Human TMIGD1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMIGD1 CRISPRa kit - CRISPR gene activation of human transmembrane and immunoglobulin domain containing 1
Symbol:	TMIGD1
Components:	GA117948G1, TMIGD1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117948G2, TMIGD1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117948G3, TMIGD1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q6UXZ0
Synonyms:	TMIGD; UNQ9372
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)