

Product datasheet for **GA117971**

Human TMEM238 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM238 CRISPRa kit - CRISPR gene activation of human transmembrane protein 238
Symbol:	TMEM238
Components:	GA117971G1, TMEM238 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA117971G2, TMEM238 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA117971G3, TMEM238 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	C9JI98
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)