

Product datasheet for **GA117983**

Human LDLRAD1 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | LDLRAD1 CRISPRa kit - CRISPR gene activation of human low density lipoprotein receptor class A domain containing 1 |
| Symbol: | LDLRAD1 |
| UniProt ID: | <u>Q5T700</u> |
| Synonyms: | low density lipoprotein receptor A domain containing 1; low density lipoprotein receptor class A domain containing 1; OTTHUMP00000010102 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)