


Product datasheet for **GA118006**

Human C22orf36 (LRRC75B) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LRRC75B CRISPRa kit - CRISPR gene activation of human leucine rich repeat containing 75B
Symbol:	LRRC75B
Components:	GA118006G1, C22orf36 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118006G2, C22orf36 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118006G3, C22orf36 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q2VPJ9</u>
Synonyms:	C22orf36; FAM211B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)