

Product datasheet for **GA118009**

Human TMEM247 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM247 CRISPRa kit - CRISPR gene activation of human transmembrane protein 247
Symbol:	TMEM247
Components:	GA118009G1, TMEM247 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118009G2, TMEM247 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118009G3, TMEM247 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A6NEH6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)