

Product datasheet for **GA118014**

Human C2orf68 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C2orf68 CRISPRa kit - CRISPR gene activation of human chromosome 2 open reading frame 68
Symbol:	C2orf68
Components:	GA118014G1, C2orf68 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118014G2, C2orf68 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118014G3, C2orf68 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q2NKK9</u>
Synonyms:	HCRCN81
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)