

Product datasheet for **GA118017**

Human PLEKHM3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLEKHM3 CRISPRa kit - CRISPR gene activation of human pleckstrin homology domain containing M3
Symbol:	PLEKHM3
Components:	GA118017G1, PLEKHM3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118017G2, PLEKHM3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118017G3, PLEKHM3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZWE6</u>
Synonyms:	DAPR; PLEKHM1L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)