


Product datasheet for **GA118018**

Human C2orf80 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C2orf80 CRISPRa kit - CRISPR gene activation of human chromosome 2 open reading frame 80
Symbol:	C2orf80
Components:	GA118018G1, C2orf80 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118018G2, C2orf80 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118018G3, C2orf80 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q0P641</u>
Synonyms:	GONDA1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)