


Product datasheet for **GA118020**

Human SNORC activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SNORC CRISPRa kit - CRISPR gene activation of human secondary ossification center associated regulator of chondrocyte maturation
Symbol:	SNORC
Components:	GA118020G1, SNORC gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118020G2, SNORC gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118020G3, SNORC gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6UX34</u>
Synonyms:	ASCL830; UNQ830
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)