

Product datasheet for **GA118030**

Human PRR23C activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PRR23C CRISPRa kit - CRISPR gene activation of human proline rich 23C
Symbol:	PRR23C
Components:	GA118030G1, PRR23C gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118030G2, PRR23C gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118030G3, PRR23C gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q6ZRP0
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)