


## Product datasheet for **GA118033**

### Human LEKR1 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                            |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| Description:  | LEKR1 CRISPRa kit - CRISPR gene activation of human leucine, glutamate and lysine rich 1                                                                                                                                                                                                                   |
| Symbol:       | LEKR1                                                                                                                                                                                                                                                                                                      |
| Components:   | GA118033G1, LEKR1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA118033G2, LEKR1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA118033G3, LEKR1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| Synonyms:     | FLJ16641; FLJ37161                                                                                                                                                                                                                                                                                         |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:



[View online »](#)