

Product datasheet for **GA118040**

Human TMPRSS11F activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMPRSS11F CRISPRa kit - CRISPR gene activation of human transmembrane serine protease 11F
Symbol:	TMPRSS11F
Components:	GA118040G1, TMPRSS11F gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118040G2, TMPRSS11F gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118040G3, TMPRSS11F gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q6ZWK6
Synonyms:	FLJ16046
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)