

Product datasheet for **GA118048**

Human PSMG4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSMG4 CRISPRa kit - CRISPR gene activation of human proteasome assembly chaperone 4
Symbol:	PSMG4
Components:	GA118048G1, PSMG4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118048G2, PSMG4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118048G3, PSMG4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q5JS54
Synonyms:	bA506K6.2; C6orf86; PAC4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)