

Product datasheet for **GA118091**

Human CSAG3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CSAG3 CRISPRa kit – CRISPR gene activation of human CSAG family member 3
Symbol:	CSAG3
Components:	GA118091G1, CSAG3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118091G2, CSAG3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118091G3, CSAG3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y5P2</u>
Synonyms:	CSAG3A; CT24.2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

