

Product datasheet for **GA118183**

Human ZNF793 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF793 CRISPRa kit - CRISPR gene activation of human zinc finger protein 793
Symbol:	ZNF793
Components:	GA118183G1, ZNF793 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118183G2, ZNF793 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118183G3, ZNF793 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZN11</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)