

Product datasheet for **GA118186**

Human PINLYP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PINLYP CRISPRa kit - CRISPR gene activation of human phospholipase A2 inhibitor and LY6/PLAUR domain containing
Symbol:	PINLYP
Components:	GA118186G1, PINLYP gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118186G2, PINLYP gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118186G3, PINLYP gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A6NC86</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)