


Product datasheet for **GA118210**

Human SPINT4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SPINT4 CRISPRa kit - CRISPR gene activation of human serine peptidase inhibitor, Kunitz type 4
Symbol:	SPINT4
Components:	GA118210G1, SPINT4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118210G2, SPINT4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118210G3, SPINT4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q6UDR6
Synonyms:	C20orf137; dj601O1.1; SPINT3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)