

## Product datasheet for **GA118211**

### Human LOC391322 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                        |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                           |
| Description:  | LOC391322 CRISPRa kit - CRISPR gene activation of human D-dopachrome tautomerase-like                                                                                                                                                                                                                                  |
| Symbol:       | LOC391322                                                                                                                                                                                                                                                                                                              |
| Components:   | GA118211G1, LOC391322 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA118211G2, LOC391322 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA118211G3, LOC391322 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                                |

#### Product images:



[View online »](#)