


Product datasheet for **GA118241**

Human TTMB (TMEM200B) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM200B CRISPRa kit - CRISPR gene activation of human transmembrane protein 200B
Symbol:	TMEM200B
Components:	GA118241G1, TTMB gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118241G2, TTMB gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118241G3, TTMB gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q69YZ2</u>
Synonyms:	TTMB
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)