


## Product datasheet for **GA118242**

### Human SLC25A35 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	SLC25A35 CRISPRa kit - CRISPR gene activation of human solute carrier family 25 member 35
Symbol:	SLC25A35
Components:	GA118242G1, SLC25A35 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118242G2, SLC25A35 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118242G3, SLC25A35 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">Q3KQZ1</a></u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)