


Product datasheet for **GA118243**

Human MEX3D activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MEX3D CRISPRa kit - CRISPR gene activation of human mex-3 RNA binding family member D
Symbol:	MEX3D
Components:	GA118243G1, MEX3D gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118243G2, MEX3D gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118243G3, MEX3D gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q86XN8</u>
Synonyms:	MEX-3D; MEX3; OK/SW-cl.4; RKHD1; RNF193; TINO
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)