

Product datasheet for **GA118271**

Human C16orf88 (KNOP1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KNOP1 CRISPRa kit - CRISPR gene activation of human lysine rich nucleolar protein 1
Symbol:	KNOP1
Components:	GA118271G1, C16orf88 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118271G2, C16orf88 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118271G3, C16orf88 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q1ED39
Synonyms:	101F10.1; C16orf88; FAM191A; TSG118
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)