


Product datasheet for **GA118307**

Human C3orf80 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C3orf80 CRISPRa kit - CRISPR gene activation of human chromosome 3 open reading frame 80
Symbol:	C3orf80
Components:	GA118307G1, C3orf80 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118307G2, C3orf80 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118307G3, C3orf80 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>F5H4A9</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)