


Product datasheet for **GA118308**

Human C4orf48 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C4orf48 CRISPRa kit - CRISPR gene activation of human chromosome 4 open reading frame 48
Symbol:	C4orf48
Components:	GA118308G1, C4orf48 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118308G2, C4orf48 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118308G3, C4orf48 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	CHR4_55
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)