

Product datasheet for **GA118320**

Human CRIP3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CRIP3 CRISPRa kit - CRISPR gene activation of human cysteine rich protein 3
Symbol:	CRIP3
UniProt ID:	Q6Q6R5
Synonyms:	bA480N24.2; TLP; TLP-A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)