


Product datasheet for **GA118322**

Human TCP10L2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TCP10L2 CRISPRa kit - CRISPR gene activation of human t-complex 10 like 2
Symbol:	TCP10L2
Components:	GA118322G1, TCP10L2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118322G2, TCP10L2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118322G3, TCP10L2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	bA517H2.3; OTTHUMP00000017617; t-complex 10-like 2; t-complex 10-like 2 (mouse)
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)