

Product datasheet for **GA118334**

Human TOMM5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TOMM5 CRISPRa kit - CRISPR gene activation of human translocase of outer mitochondrial membrane 5
Symbol:	TOMM5
Components:	GA118334G1, TOMM5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118334G2, TOMM5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118334G3, TOMM5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8N4H5
Synonyms:	bA613M10.3; C9orf105; Tom5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)